



JADA 2.0 USER GUIDE

User manual based on the JADA 2.0 prototype

-  <https://www.dauvea.it/>
-  <https://it.linkedin.com/company/dauvea>
-  https://twitter.com/dauvea_it
-  https://www.instagram.com/dauvea_/
-  <https://www.facebook.com/dauveasrl>



Document information

Document ID:	
Status	<input type="checkbox"/> DRAFT <input type="checkbox"/> REVISED <input checked="" type="checkbox"/> APPROVED
Classification level	<input checked="" type="checkbox"/> PUBLIC <input type="checkbox"/> PRIVATE <input type="checkbox"/> CONFIDENTIAL

Version history

Version	Date	Written by	Status
1.0	29/01/2025	Beatrice Porcu	DOC DRAFT, SECTIONS 1-2
2.0	30/01/2025	Beatrice Porcu	SECTIONS 3-5
3.0	31/01/2025	Beatrice Porcu	SECTION 6
4.0	02/03/2026	Davide Rubiu	SECTION 5

Revision history

Reviewer	Approved by	Date	Version	Review notes
Roberta Lecca Giacomo Anedda	Roberta Lecca Giacomo Anedda Beatrice Porcu	31/01/2025	1.0_FINAL	Review complete



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
 Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
 Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



INDEX

Document information	1
Version history	1
Revision history	1
JADA VERSION 2.0 – USER GUIDE	3
1. ACTIVATING YOUR ACCOUNT	4
2. LOGGING IN	6
3. HOME PAGE	7
a. ADMIN AND COACH ONLY	7
b. ATHLETES	7
4. WORKING WITH YOUR SENSORS (ADMIN AND COACH ONLY)	8
a. SAVE A NEW SENSOR	8
b. CONNECT AND DISCONNECT AN EXISTING SENSOR	9
c. DELETE A SENSOR	10
5. SESSIONS	11
a. SESSION DASHBOARD	11
b. REGISTER A NEW SESSION (ADMIN AND COACH ONLY)	12
c. YOUR SESSION'S RESULTS (ADMIN AND COACH ONLY)	15
d. SEARCH IN YOUR SESSION HISTORY	15
e. DELETING A SESSION (ADMIN AND COACH ONLY)	17
6. YOUR SETTINGS	18
a. YOUR DETAILS	19
b. YOUR SECURITY DETAILS	19
c. YOUR TEAM PAGE	20
d. CREATE A NEW USER (COACH AND ADMIN ONLY)	21
e. MODIFYING AN EXISTING USER (COACH AND ADMIN ONLY)	23
f. YOUR TRAINING GROUPS	24



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



JADA VERSION 2.0 – USER GUIDE

The guide explains the principal user flows for the first version of JADA. All information and screenshots are to be considered valid for all user types (admin – coach – athlete) unless otherwise specified.



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it




1. ACTIVATING YOUR ACCOUNT

- 2. In the inbox of your email account, you'll find an email from 'supporto@dauvea.it' with the object "Activate your JADA account".

To: <user@guide.com>

[Activate your JADA account](#)

- 3. In the email, click the link to activate your account.



- 4. In the verification page, you can choose to register your account via email and a custom password that you choose, or verify your account via Microsoft SSO or Google.


Welcome to Jada


You can access the app by creating a personalized password for your account or via the SSO of your associated email.

Log in with email and password

EMAIL & PASSWORD

Log in via SSO

 Sign in with Google


 **MICROSOFT**

BEWARE! To insert your custom password you can click the button "EMAIL & PASSWORD" to show the registration form.

Log in with email and password

EMAIL & PASSWORD

user@guide.com

Password 

Confirm Password

Passwords must match

LOGIN



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



5. If you decide to verify your account via custom login, you will be required to add a second authentication via a 2FA app (e.g. Microsoft Authenticator or Google Authenticator).

Activate 2FA

For your account safety, we require you to link your account to a 2FA app. We recommend either Google Authenticator or Microsoft Authenticator

STEP 1: Scan code with your app



STEP 2: Insert below the code produced by your authenticator app

Can't scan the QR code? Copy this string in your app

6QSI32EFNPQP7AZ4

6. Once the verification process is completed with successful completion, you will be redirected to the login page to log into the application.




Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it




2. LOGGING IN



1. Log in at <https://jada-project.dauvea.it/> using the method you chose during the verification process.

2. If you log in via custom login, you will be required to insert the OTP code from your authentication app.




Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



3. HOME PAGE


a. ADMIN AND COACH ONLY



In the Home page you can find a synthetic review of all you app sections.

1. **ATHLETES SECTION:** here you can find a comprehensive list of all your athletes, either active or not. Clicking on the arrow next to the label will take you to the Team's settings (for details, see the dedicated section)
2. **SENSORS SECTION:** here you can see a list of the previously registered sensors. Clicking on the arrow next to the label will take you to the Sensors page (for details, see the dedicated section)
3. **SESSION SECTION:** the graph shows all sessions registered in the last 30 days. If you don't have registered any session yet, the graph will not be shown. As with the previous sections, clicking on the arrow next to the label will take you to the Session dashboard (for details, see the dedicated section)
4. **MENU BUTTON:** clicking on the hamburger button on the left of the header will open the side menu where you will find links to all key features of the application, a logout function and a link to the Changelog page, where you will find more information on the latest release of the JADA app.

b. ATHLETES



The athletes' home page is slightly different. It contains only the data pertaining to the athlete who logged in (team info and session registered by them)




Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it




4. WORKING WITH YOUR SENSORS (ADMIN AND COACH ONLY)

This is your sensors' home. From this page you can save new sensors, connect your existing ones or delete them.



Before you start working with your devices, you need to activate some features in your browser. The process is quite simple and you will find a tutorial similar to this one from inside the application clicking on the HELP button on top of the page:

- From your Chrome browser, go to the page `chrome://flags/`
- In the page that opens, search and activate the following flags:
 - o `#enable-experimental-web-platform-features`
 - o `#enable-web-bluetooth-new-permissions-backend`: set it either to default



a. SAVE A NEW SENSOR


1. Click the button SAVE NEW SENSOR. A modal will appear where you can choose to go directly to the connection process or to see a brief tutorial that will guide you through the process.




Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it





2. Once the connection is established, you will see the sensor data and you can save the device in your sensor list by clicking the SAVE button.



3. Once the device is saved and renamed based on your team initials, you will see it connected in your list. You can identify the connected devices by the orange Bluetooth icon and the battery value indicator. When the sensor is disconnected instead, the icon is dark blue and no battery value is shown.



b. CONNECT AND DISCONNECT AN EXISTING SENSOR

You can connect and disconnect a previously saved device by clicking the Bluetooth icon shown above.

On mobile the Bluetooth button is a normal button with the same icon and underneath the sensor data, but the functioning is exactly the same.



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



c. DELETE A SENSOR

You can delete a sensor by clicking on the trash can button. If you confirm the deletion, the sensor will be renamed to the factory name and deleted from your sensor list. If you want to add it again, you can follow the instructions to save a new sensor again.




Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



5. SESSIONS

a. SESSION DASHBOARD



- Sections 1 and 3: some statistics on the sessions you registered in the past. On 1 you can see your sessions divided by type, on 3 you can see the athletes with most sessions registered
- Section 2: last 5 sessions saved; to see and search among your other sessions, click on 5
- Section 4: START NEW SESSION button; click to start a new session
- Section 5: SESSION HISTORY button, where you can search among all your registered sessions



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
 Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
 Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



b. REGISTER A NEW SESSION (ADMIN AND COACH ONLY)

The new session wizard is accessible via the session dashboard described above, or via the sidebar menu by clicking 'New session'.

The wizard will guide you through the process.

1. Select your sensor. If it's disconnected, the browser's popup will allow you to connect it. Once your sensor has blinked in yellow and has made a sound, you can go forward.

Step 1/3. Select sensor to use for the session

Sensors

<input checked="" type="radio"/> TM2-4C65 - 221e	<input type="radio"/> DVL-4B08 - 221e
<input type="radio"/> DVL-4254 - 221e	<input type="radio"/> PLL-79FE - 221e
<input type="radio"/> PLL-AC5C - 221e	<input type="radio"/> PLL-C567 - 221e
<input type="radio"/> PLL-AA12 - 221e	<input type="radio"/> PLL-F6BF - 221e
<input type="radio"/> PLL-5235 - 221e	<input type="radio"/> PLL-8C1C - 221e
<input type="radio"/> PLL-52F4 - 221e	<input type="radio"/> PLL-6A81 - 221e
<input type="radio"/> PLL-8DEF - 221e	<input type="radio"/> PLL-793F - 221e
<input type="radio"/> PLL-8454 - 221e	<input type="radio"/> PLL-8027 - 221e
<input type="radio"/> PLL-96D8 - 221e	<input type="radio"/> PLL-7B24 - 221e

FORWARD

2. Select your athlete. You can select your athlete only among those that activated their account.

Step 2/3. Select athlete to monitor

Athletes

<input type="radio"/> Alessandra Concu	<input type="radio"/> Marta Pusceddu	<input type="radio"/> Alessia Panzeri
<input type="radio"/> Valentina Biggio	<input type="radio"/> Rossella Pittalis	<input type="radio"/> Irene Muru
<input type="radio"/> Viola Zinzula	<input type="radio"/> Beatrice Carta	<input type="radio"/> Tecla Puddu
<input type="radio"/> Ginevra Murgia	<input type="radio"/> Asia Melis	<input type="radio"/> Atleta 12

BACK FORWARD




Dauvea s.r.l.


Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



3. Select session type: the available sessions are Posturography (static), Jumps or Time To Stability. Training sessions are coming in a future release. The forward button will activate once you select your type.



4. Select sub session type, based on the session type chosen before. The forward button will activate once you select your sub type.



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



5. In case you select a personalized duration, you will have to compile the time duration forms that appears underneath the list, and click the save button to save the personalized duration.

Step 1 Step 2 Step 3 **Step 4**

Step 4/4. Select session type

Available session types:

Session type: Monopodalic Bipodalic

Eyes: Eyes closed Eyes open

Duration: Predefined duration- 30 sec. Custom duration

Custom duration (mm:ss)

0 : 30

BACK START NEW SESSION



Once you receive the correct saving message you can start a new session.

Loading session



00:03:60




Dauvea s.r.l.
Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



c. YOUR SESSION'S RESULTS (ADMIN AND COACH ONLY)

At the end of the countdown, you will see a recap of the session you just registered.



The recap contains:

- Your athlete's data
- The session's info (type, frequency, duration)
- Your action buttons:
 - o DELETE deletes the session info and all sensor data registered during it
 - o DOWNLOAD exports a file with the raw sensor data, useful for debugging purpose
 - o SAVE saves your session and adds it to your session history to be accessed at a later date
- Benchmark drawer: opens to show your session data compared to a gold standard based on literature

Benchmark

Posturography - monopodalic posturography - eyes open - personalized duration

swayComputeIMU - ellipseArea: 0.0036



swayComputeIMU - totalPathLength: 1.4643



swayCompute - ellipseArea: 41.2461



swayCompute - totalPathLength: 159.8438




- Session graphs drawer: opens to show your session graphs

Session completed



Once you confirm or delete your session, you can start over by registering another one or going back to your dashboard.



d. SEARCH IN YOUR SESSION HISTORY

You can access your sessions' history by the button on the session dashboard or from the side menu. The page presents at the top the available filters you can choose



Dauvea s.r.l.



Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



from by clicking on each one of them to expand it, a button to start your search and the result of your search at the bottom.

You can filter by:

- Date interval (if an end date is not specified, it filters to today's date)
- Your team's athletes
- Session types
- Coaches who registered the session



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



e. DELETING A SESSION (ADMIN AND COACH ONLY)

Other than during the registration process, you can also delete a session from your session history clicking on the three dots on the right of the session card and clicking on the DELETE button.

Posturography - Monopodalic posturography
- eyes open - personalized duration

Coach Test

Giulia Bianchi

100Hz - 5 sec.



1/30/25

swayComputeIMU - ellipseArea: **0.0036** (-0.5148 - 0.9770)

swayComputeIMU - totalPathLength: **1.4643**(1.3352 - 46.3986)

swayCompute - ellipseArea: **41.2461** (-4,938.5525 - 9,400.6615)

swayCompute - totalPathLength: **159.8438** (20.7374 - 4,682.2636)




Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it




6. YOUR SETTINGS

You can access your team data either from the Home page from the button on top of your athletes list (for coaches and admins only) or from the Settings in the sidebar menu that you access from the menu button on the left of your header bar.



This is your settings area. From here you can access your account's data and a bird's eye view on all team members.

You can move between the different features by clicking on the tabs underneath your email address.



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



a. YOUR DETAILS

Personal details

Coach Test

Preferred language

IT
 EN

SEND

Here you can change your personal data. Athletes can also set their height and weight in this page.

b. YOUR SECURITY DETAILS

Authentication settings

Authentication type: Custom

MODIFY RESET PASSWORD

DELETE ACCOUNT

Here you can see more information on your authentication method, reset your password if you are

logging in with a custom form, or delete your account. You can also change your authentication method by clicking the MODIFY button. If you choose to do so and change the authentication method to a custom one, you will receive an email notification to set your password.

Set auth type

You can change your preferred type of authentication to log in, **but your email address will remain the same.**

MICROSOFT GOOGLE CUSTOM



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



c. YOUR TEAM PAGE

The team page contains all important info and statistics pertaining your team. By opening and closing the drawers underneath the team's statistics, you can see the list of team members divided by role.

CT Coach Test
sviluppo@dauvea.it

DETAILS LOGIN & SECURITY **TEAM** TRAINING GROUP

Team 2
www.palestra1.it
emailpalestra1.test.it
0704444444

Active users: **48%**

Admin: **1** Coaches: **4** Athletes: **26**

REGISTER NEW USER

Admin ▾
Coach ▾
Athletes ▾



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



d. CREATE A NEW USER (COACH AND ADMIN ONLY)

1. From the team details section in your settings page, click on the “Register new user” button.

The screenshot displays the 'TEAM' tab for a user named 'Coach Test' (sviluppo@dauvea.it). The interface includes a navigation bar with tabs for 'DETAILS', 'LOGIN & SECURITY', 'TEAM', and 'TRAINING GROUP'. The 'TEAM' tab is active, showing 'Team 2' with contact information: www.palestra1.it, emailpalestra1.test.it, and 0704444444. A summary card indicates 48% active users, 1 Admin, 4 Coaches, and 26 Athletes. A 'REGISTER NEW USER' button is visible, along with a dropdown menu for user roles: Admin, Coach, and Athletes. A large 'REGISTER NEW USER' button is highlighted at the bottom of the screenshot.



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



- In the form that opens, insert the new user data. Keep in mind that the user will need a valid and unique email address to receive the email to activate their account.
ONLY FOR ADMINS: you can also register coaches using the same form and selecting “coach” as user role.

Register new user X

Personal details

Name Last name

Email

Height in meters Weight in kg

Select gender

Jan 30, 2025 Birth date

Preferred language

IT



- Once all data is correctly inserted, click on the SEND button.




e. MODIFYING AN EXISTING USER (COACH AND ADMIN ONLY)

1. You can access the form to modify a user from the Home page or the team settings page, clicking on the three dots next to the user name. From there you can see the user details or modify the user's data.

Giulia Bianchi - giulia.bianchi@gmail.com



2. Click on the MODIFY button. A form with the user's info will open. Modify all the necessary fields, then click SEND at the bottom of the form.



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it



f. YOUR TRAINING GROUPS

Here you can access all training groups you created, delete them or create new ones.

Training groups will be available to use in future implementations, but you can already create them now.

Team's training groups

The screenshot shows a grid of six training group cards. Each card contains a red circular icon with a white plus sign, followed by the group name and its creation date: trainingTest_06 - 11/11/24, trainingTest_05 - 11/6/24, trainingTest_04 - 11/6/24, trainingTest_03 - 11/6/24, trainingTest_02 - 11/6/24, and trainingTest_01 - 11/6/24. Below the grid is a prominent red button labeled 'CREATE NEW GROUP'.

The form is titled 'Training group name' and features a text input field. Below this is a section titled 'Athletes' with a list of names, each accompanied by an unchecked checkbox and a small circular icon with a plus sign. The names listed are: Giulia Bianchi, Chiara Ferrari, Elisa Russo, Valentina Marini, Sofia Conti, Sara De Luca, Martina Verdi, Francesca Romano, Alessia Galli, Beatrice Moretti, Ilaria Guerra, and Testmodifica rweASDA. A vertical red scrollbar is visible on the right side of the athlete list.

To create a new group, click on the CREATE NEW GROUP BUTTON and in the form that appears insert a name (mandatory), select your athletes and choose if you want the other team members to see your training group (yes by default).

Make group visible to the rest of the team

SAVE

Once you click SAVE the training group will be available to use.



Dauvea s.r.l.

Capitale sociale 10.000€ i.v. - P.IVA 03718890928
Sede operativa: via Cesare Battisti, 14 - 09123 Cagliari - Italy
Sede legale: via Alberto Riva Villasanta, 233 - 09134 Cagliari - Italy
info@dauvea.it - www.dauvea.it

